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Importance of Rheology and Viscosity in Paper and Board Coating

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Abstract

Market demand for coated graphical paper has declined steady as a result of electronic media. However, other coated grades such as coated board and barrier coated packaging papers and specialties are forecasted with decent growth grades. Despite of the actual market situation the operational environment and coating functionalities are under continuous investigation. Flow properties i.e. rheology plays an important overall role in optimizing runnability and quality. This paper demonstrates the importance of rheology and viscosity in paper and board coating. A focus is in high shear viscosity that has been demonstrated to be in a critical parameter in coater runnability. Its importance will be discussed both from theoretical and practical standpoints. Furthermore, a new instrument for measuring high shear viscosity fast and easy will be discussed.

Keywords : High shear measurement, Rheological properties, Viscosity, Capillary viscometer, Coating color, Coater runnability

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Suimeidai, Kawanishi-shi, Hyogo 666-0116, Japan

^{*2}E-mail : info@snc94.co.jp


Fig. 1 Structural and non-structural viscosity.

Table 1 Governing factors in low and high shear viscosity

Surface chemistry dominates in low shear viscosity	Hydrodynamics dominate in high shear viscosity
- Van der Waals - Electrostatic Repulsion/Attraction - Steric Factors	- Particle size - Distribution of Particle size - Shape of Particles - Viscosity of Water Phase - Solid Fraction
Can be measured with Brookfield, Anton Paar and Hercules rheometers	Can be measured with capillary viscometer


Fig. 2 Shear rate in viscosity measurement to simulate the real coating process condition.

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
Fig. 3 Effect of filler particle size distribution (PSD) on viscosity (measured by ACA viscometer).

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- Same average particle size.
- Narrow PSD gives rise to higher viscosity.
- More blade streaks and bleeding with narrow PSD pigment.
- Uneven blade loads with narrow PSD pigment.


Fig. 4 Industrial example of ideal viscosity and deviation.



Fig. 5 Viscosity and blade load (runnability). Non-structural viscosity has a very high correlation with the actual coating process.

Brookfield (low shear): Viscosity is the same regardless of the size of latex particles.

Capillary (high shear): Viscosity decreases strongly with decreasing particle size.

- Smaller particles fit easier between pigment particles
- Possible to increase solids with smaller particle size latex

Fig. 6 Effect of latex particle size on viscosity.

Table 2 Effect of thickeners on viscosity

Property	Viscosity(mPa·s)		
	Without thickeners	CMC	Starch
Viscosity of water phase	1	3.5	8.3
Brookfield viscosity	100	1200	750
High shear viscosity	12	40	55

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Reference

- 1) Esa Lehtinen : Papermaking Science and Technology, Book 11, Pigment Coating and Surface Sizing of Paper, 661, Gummerus Printing (2000)